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(TI-010)

Dental Zirconia

HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered)

ENGLISH INSTRUCTIONS FOR USE


I. Introduction

KATANA Zirconia has 2 shapes available: Block and Disc. This IFU is for KATANA Zirconia HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered). KATANA Zirconia HT & ML /HTML are a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia HT has 6 available thicknesses: (10mm, 14mm, 18mm, 22mm, 26mm and 30mm) and has 3 shade variations: (HT10, HT12 and HT13). KATANA Zirconia ML has 3 available thicknesses: (14mm, 18mm and 22mm) and has 6 shade variations: (A Light, A Dark, B Light, C Light, D Light and A White). KATANA Zirconia HTML has 3 available thicknesses: (14mm, 18mm and 22mm) and has 14 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia HT is higher in translucency and KATANA Zirconia ML / HTML consists of 4 graded shade layers. KATANA Zirconia HT & ML / HTML are recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).


III. Sintering Program

Temperature	Programming Rate	Holding Time
Room Temp. — 1500°C (2732°F)	10°C/min. (-18°F/min.)	—
1500°C (2732°F)	—	2 hrs
1500°C (2732°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

Sintering Program 2

Temperature	Programming Rate	Holding Time
Room Temp. — 1515°C (2759°F)	35°C/min. (63°F/min.)	—
1515°C (2759°F)	—	0,5 hrs
1515°C (2759°F) — Room Temp. (*)	-45°C/min. (-81°F/min.)	—

* The restorations may be removed at 800°C (1472°F) or less depending upon circumstances.

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

- There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2 is compatible.
- Sintering is only recommended for crowns, veneers, and inlays/onlays with Sintering program 2.
- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.
- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercut, guide grooves, the formation of retentive holes, and sharp corners.
- Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Anterior crown or bridge	0.4 mm or more
Veneer	0.4 mm or more
Posterior crown or bridge	0.5 mm or more
Inlay or onlay	0.5 mm or more

12. Use the following cross-sectional areas for connectors when fabricating bridges.

Location & indication	Connector cross section
Anterior	2- or 3-unit bridges 7 mm ² or more <p>more than 4-unit bridges 9 mm² or more</p>
Posterior	2- or 3-unit bridges 9 mm ² or more <p>more than 4-unit bridges 9 mm² or more</p>

- Keep the pontic in the bridge within 2 teeth. When the pontic of 2 teeth continues, keep the connector cross section between the pontic at 12 mm².
- Keep the cantilever bridge to 1 pontic tooth and the connector cross section at 12 mm².
- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations. When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.
- Caution in conjunction with usage:**
 - Do not use the disc if there is a crack noticed after removing it from the package.
 - Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
 - If you find a crack in the restorations, do not use.
 - The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
 - When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
 - Cut and remove the plastic ring prior to sintering the whole disc.

- Do not use the disc if there is a crack noticed after removing it from the package.
- Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
- If you find a crack in the restorations, do not use.
- The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
- When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
- Cut and remove the plastic ring prior to sintering the whole disc.

Storage:

- Store in a cool and dry place. Keep away from direct sunlight.
- The product should NOT be stored at 10-30°C (50-86°F).
- Do not remove the disc from its packaging during storage.
- The disc is fragile, and requires care when handling.
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

Caution:

- This product should NOT be used when malocclusion, clenching or bruxism conditions are present.
- When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.

Warning:

- If the patient is hypersensitive to zirconia or any other components, this product not be used.

Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

Contraindications:

- If the patient is hypersensitive to zirconia or any other components, this product not be used.

Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

Caution:

- This product should NOT be used when malocclusion, clenching or bruxism conditions are present.
- When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.

Caution:

- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated by the furnace with your bare hands.
- Do not take the restorations out of the sintering furnace during high temperature, as te quenching causes the breaking. However, if an elevating type furnace is used under the sintering program 2, the restorations may be removed from the furnace at 800°C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.

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